


Introduction on Patient Monitoring and Wearable Device Reference Designs

Agenda

- TI in medical
- Market trend
- Fundamentals and challenges
- Reference Designs (TIDA-01614, TIDA-010005 and etc)

TI Semiconductors in every medical category

Medical Imaging



Patient Monitoring & Diagnostics



Medical Equipment



Home Healthcare



Personal Care & Fitness





Patient Monitoring Market Trend


- The growing elderly population drives the high demand of remote monitoring- In EU, 700B Euro spent on chronic disease per year
- Remote monitoring enhances quality of care and reduces healthcare cost
- Wearable wireless medical technology enables accurate and reliable data in a smaller form factor: multi-modalities, longer battery, SHIP mode
- Artificial Intelligence uses analytics and big data to improve decision making and early prevention

Patient Monitoring Basics

The electrocardiogram (ECG)
measure of electrical activity of the heart




Photoplethysmography (PPG) is an *optical* measurement of an organ's volume.



Patient Monitoring Basics

- **Pulse Transit Time (PTT)** – time interval between the R-peak of the ECG and the max slope of the PPG
- Involves simultaneous ECG and PPG measurements
- Systolic blood pressure (SBP) can be estimated from PTT



ECG vs. PPG


ECG	Feature Description	PPG
Electrical	Measurement type	Optical
Electrodes	Sensor type	Photodiode
Yes	Can measure heart rate?	Yes
Yes	Diagnostic Information	Yes
2 (Across chest)	Minimum number of skin contacts required?	1 (Finger or wrist)
≥ 1	Number of ADC channels required	1

ECG Lead and ADC Channels

IEC60601-2-51 – Diagnostics


Number of Leads	Leads Used	Number of ADC Channels
1	Lead I	1
3	Lead I, Lead II, Lead III	2
6	Lead I, Lead II, Lead III, aVR, aVL, aVF	2
12	Lead I, Lead II, Lead III, aVR, aVL, aVF, V1 – V6	8

Standards	Electrodes Needed
1 Lead	LA, RA
3 Lead	LA, RA, LL
6 Leads	LA, RA, LL
12 Leads	LA, RA, LL, V1-6




ECG characteristics

Frequency domain



Challenges in measuring ECG



Challenges in Optical Bio-Sensing

- Low power for longer battery life



- Skin tone variation



- Best PPG signal for Motion cancellation Algorithms



- Performance with Glass

- Low temperature performance



- Ambient Light



Medical Sector Page on TI web

The screenshot shows the Texas Instruments (TI) website's medical sector page. The top navigation bar includes links for Products, Applications & designs (which is highlighted with a blue border), Tools & software, Support & training, Order Now, About TI, My History, Cart, English, and myTI. The main content area is organized into several sections:

- Imaging:** CT and PET scanners, MRI, Ultrasound scanners, X-ray systems.
- Patient monitoring and diagnostics:** Clinical digital thermometer, Digital stethoscope, Electrocardiogram (ECG), Endoscopes, Eye, ear, nose and throat exam, In vitro diagnostics, Medical sensor patches, Mother and neonatal care monitors, Multiparameter patient monitor, Pulse oximeter, Sleep diagnostics.
- Home healthcare:** Blood glucose monitors, Blood pressure monitors, CPAP machines, Electronic thermometers, Hearing aids, Oxygen concentrator, Telehealth systems.
- Medical equipment:** Anesthesia delivery system, Dialysis machines, Electronic hospital bed & bed control, Infusion pumps, Medical accessories, Medical chairs and tables, Motorized electronic wheel chair, Surgical equipment.
- Personal care and fitness:** Electric toothbrush, Fitness machines, Wearable fitness & activity monitor, Beauty and grooming.
- Central Content Area:** Displays a thumbnail image of a circuit board and the URL [medical/overview.html](#). Below the image are links for Corporate systems >, Industrial >, Personal electronics >, and Reference designs >. The **Medical** category under Industrial is highlighted with a blue border and an arrow points to it from the right sidebar.
- Right Sidebar:** Includes links for Aerospace & defense, Building automation, Electronic point of sale (EPOS), Factory automation & control, Grid infrastructure, Industrial transport (non-car & non-light truck), Lighting, Motor drives, Power delivery, Pro audio, video & signage, Test & measurement, Connected peripherals & printers, Data storage, Gaming, Home theater & entertainment, Mobile phones, PC & notebooks, Portable electronics, Tablets, TV, and Wearables (non-medical). It also features a note about reference designs and an Internet of Things link.

Medical Sensor Patches

Temperature Analog Front End

Temperature AFE drives the temperature sensors as well as conditions the signal coming from analog temperature sensors. The temperature signal is converted to digital domain by using ADC for further processing.

[View more](#)

REFERENCE DESIGNS (11)

PRODUCTS (15)

Sensors (4)

Analog Temperature Sensors (1)

Digital Temperature Sensors (3)

TMP117 – $\pm 0.1^\circ\text{C}$ accurate digital temperature sensor with integrated NV memory

TMP112 – 1.4V-Capable $\pm 0.5^\circ\text{C}$ Accuracy Digital Temperature Sensor in the Compact SOT-563 Package

TMP102 – 1.4V-Capable Temperature Sensor with I²C/SMBus Interface and Alert Function in SOT-563

Data converters (5)

Switches & multiplexers (2)

Amplifiers (2)

Power management (2)

port & training

- Simplify integration
- Achieving low power
- Processing

Wi-Fi

Sub GHz

RF Interface

Memory Card

AMP

Speaker

User Interface

PPG/Optical/Spectroscopy Front End

The optical front end has two sections: one for driving the LEDs and other for processing the signal received from photo diodes. LED drivers are operated with digital data coming from DAC and transimpedance amplifier with a PGA conditions signal coming from photo diode. This subsystem can be used for PPG as well as spectroscopy measurement.

[View more](#)

UCTS (31)

t monitor

REFERENCE DESIGNS (7)

PRODUCTS (31)

Amplifiers (8)

(voltage (11) / current (4))

General-Purpose Op Amps (3)

Instrumentation Amplifiers (1)

Switches & multiplexers (6)

Data converters (14)

Biosensing AFEs (6)

AFE4420 – Ultra-small integrated AFE with FIFO for multisensor wearable optical heart-rate monitoring

AFE4900 – Ultra-low-power integrated AFE for wearable optical, electrical biosensing

AFE4403 – Ultra-Small Integrated Analog Front End (AFE) for Heart-Rate Monitors and Low-Cost Pulse Oximeters

AFE4400 – Integrated Analog Front End (AFE) for Heart-Rate Monitors and Low-Cost Pulse Oximeters

AFE4490 – Integrated Analog Front End (AFE) for Pulse Oximeters

Multiparameter Patient Monitor

Multiparameter patient monitor

Multiparameter Patient Monitor integrated circuits and reference designs

Description | Reference designs & products | Technical documents | Support & training


Description

Our integrated circuits and reference designs help you create highly accurate multiparameter patient monitors to meet system design requirements while minimizing power consumption.

Innovative multiparameter patient monitors require:

- Precise, multi-parameter measurements with high sensitivity
- Wireless, secure connectivity for data logging and monitoring
- Intuitive, responsive human machine interface (HMI)

Select a subsystem for: Patient monitor - ECG module



Analog front end

The AFE digitizes the heart signal after signal conditioning the signal. It includes multiplexer, instrumentation amplifier, filter, ADC. Additionally, it includes temperature sensor, respiratory circuit, pacemaker detection and lead off detection circuit.

[View more](#)

REFERENCE DESIGNS (9)

PRODUCTS (50)

Multiparameter front end reference design for vital signs patient monitor



- [Schematic/Block diagram](#)
- [Reference guide](#)
- [View reference design](#)

Alarm tone generator reference design

Software-Configurable Cardiac Pacemaker Detection Module Reference Design



Reference Design

Design Guide: TIDA-01614
Multiparameter Front-End Reference Design for Vital Signs Patient Monitor



Description


This reference design is for a multiparameter front-end of a patient monitor that measures vital sign parameters like electrocardiogram (ECG), heart rate, SpO₂, and temperature. It uses discrete front-end integrated circuits, like the AFE4403 and ADS1292R devices, to measure these parameters. The design also uses three TMP117 sensors to accurately measure skin temperature. The design can interface with the pace detection module to detect the pace pulse. The design also includes an RS-232C/UART connection to transfer data to a computer. The entire front-end subsystem runs on a rechargeable 3.7-V Lithium-ion (Li-Ion) battery.

Resources

TIDA-0161
AFE4403
BQ24232
MSP430F4
ADS1292R
TMP117
TIDA-0100

Comprehensive Design Guides

- Pulse oximeter
- Electrocardiogram (ECG)



Boards Available for Evaluation

The following list provides details about the design:

- Supply Voltage = 5 V
- Charging current = 0.1 A

Detailed Design Considerations & Applications Info

- Monitors ECG, heart rate, SpO₂ %, respiration rate, and skin temperature
- Uses bio-sensing front-end AFE4403 for SpO₂ and heart rate measurement and ADS1292R for ECG and respiration measurement
- Includes three red LEDs and three photodiodes with ambient subtraction to improve signal-to-noise Ratio (SNR) for SpO₂ and heart measurement
- Single lead ECG Measurement with RLD
- Survives three 0.1 Cal/sec annular sensors to

$$RISET = R99 = [870 \text{ }\Omega / 0.1 \text{ A}] = 8.7 \text{ k}\Omega$$

How the input current limit (ILIM) is set:

- RLIM = KILIM / I1-MAX
- KILIM = 1530 $\text{A}\Omega$

Bill Of Materials

Designator	Quantity	Value	PartNumber	Manufacturer	Description	PackageReference
IOP31	1	0.1uf	TIDA-01614	Any	Printed Circuit Board	0402
C1, C2, C4, C6	12	0.1uf	C1005XTR1E104M050BB	TDK	CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
C3, C5, C7, C8, C9, C10, C11	1	0.01uf	C1005XTR1E104M050BB	TDK	CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C2, C3, C5, C7	4	0.01uf	GRM		CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C4	1	0.01uf	C1005XTR1E104M050BB	TDK	CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C5	1	0.01uf	C1005XTR1E104M050BB	TDK	CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C7, C9, C10	2	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
C11, C13, C24	6	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
C25, C26, C27	2	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
C15, C17	2	0.01uf	GRM		CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C16	1	1500uf	GRM		CAP, CERAM, 1500 μF , 25%, X7R, 0402	0402
C17, C19, C20, C21, C22, C23	5	0.1uf	0603		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
TR1	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
TR2	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
C18, C20	2	2200uf	C105		CAP, CERAM, 2200 μF , 25%, X7R, 0402	0402
C21	1	0.01uf	C1005XTR1E104M050BB	TDK	CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C13, C17, C22	2	2200uf	C105		CAP, CERAM, 2200 μF , 25%, X7R, 0402	0402
C23	1	0.01uf	C1005XTR1E104M050BB	TDK	CAP, CERAM, 0.01 μF , 25%, X7R, 0402	0402
C27, C30	2	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
SR_0003	1	1uf	GRM		CAP, CERAM, 1 μF , 25%, X7R, 0402	0402
ENP1	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP2	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP3	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP4	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP5	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP6	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP7	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP8	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP9	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP10	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP11	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP12	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP13	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP14	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP15	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP16	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP17	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP18	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP19	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP20	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP21	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP22	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP23	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP24	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP25	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP26	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP27	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP28	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP29	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP30	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP31	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP32	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP33	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP34	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP35	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP36	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP37	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP38	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP39	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP40	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP41	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP42	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP43	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP44	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP45	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP46	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP47	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP48	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP49	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP50	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP51	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP52	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP53	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP54	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP55	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP56	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP57	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP58	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP59	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP60	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP61	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP62	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP63	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP64	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP65	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP66	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP67	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP68	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP69	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP70	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP71	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP72	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP73	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP74	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP75	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP76	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP77	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP78	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP79	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP80	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP81	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP82	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP83	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP84	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP85	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP86	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP87	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP88	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP89	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP90	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP91	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP92	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP93	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP94	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP95	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP96	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP97	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP98	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP99	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP100	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP101	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP102	1	0.1uf	GRM		CAP, CERAM, 0.1 μF , 25%, X7R, 0402	0402
ENP103	1	0.1uf	GRM		CAP, CERAM, 0.1 $\mu\text{F$	

Features


- System measures ECG, Heart Rate, SPO₂, Respiration rate using ADS1292R and AFE4403 and Skin temperature using TMP117
- Circuit enables three electrode operation including right leg drive with good CMRR
- Pace detection circuit indicates presence of pacemaker
- Supports three 0.1 Celsius accurate sensors (TMP117) to measure the skin temperature
- Enables data transfer over isolated UART interface
- Works with 3.7V Li ion rechargeable battery
- On board memory for data logging


Applications

- Multiparameter Patient Monitor
- Medical Sensor Patches
- Pulse Oximeter
- Electrocardiogram (ECG)

Benefits

- Single IC does both ECG ,Respiration.
- Pace Detection
- ECG with 3 electrodes
- Three temperature sensors for temperature measurement





Design Challenges TIDA-01614 Solves

Design Challenge 1:


Integration of multiple modalities at optimum SNR levels and small form factor

- Monitoring of ECG, Heart Rate, SPO₂, PTT, Respiration rate and Skin temperature
- Single Lead ECG with RLD (ADS1292R)
 - > Signal amplitude: 0.2mV~2mV (p-p);
 - > BW 0.05 Hz to 2000 Hz
- Supports 3 LED and 3 Photodiodes with ambient subtraction for SPO₂ and Heart Rate monitoring with AFE4403
- Supports three 0.1 Celsius accurate sensors to measure the skin temperature (TMP117)


Design Challenge 2:

Protection and isolation against ESD & defibrillation per IEC 60601-1-2

- Using TVS0500 and series current limit resistor. TVS0500 has lower clamping voltage to protect ADS1292R.



- Isolated UART interface using an onboard MSP432P401, ISOW7842, TRS3232






ECG Analog Front End

TIDA-01614: Multiparameter front-end reference design for vital signs patient monitor



Important Parameters:


- Input Bias Current
- Input Impedance
- Input Current Noise
- Input Voltage Noise
- Power Consumption
- DC/AC CMRR



TIDA-01614 Test Setup and Test Results

Design Specs

CHARACTERISTICS	SPECIFICATIONS
ECG	One lead ECG operation with RLD. Sampling rate of 500 samples per second, supports ECG sensitivity of 100 μ V
SPO2 Measurement	Works in transmissive SPO2, refresh rate of 500 Hz
Skin Temperature Measurement	Three temperature sensor with 0.1 degree accuracy
Pace pulse Rise-time (TR) measurement range	30–200 μ s
Pace pulse duration (TD) measurement range	0.1–2 ms
Input Pace signal amplitude range	8 mV–700 mV
Input Voltage (Vin)	5 V from Micro-USB



Test Setup




GUI Display


ADS1291/2

ADC specifications

- Low-noise, high input impedance front end PGAs
- 24-bit simultaneous sampling delta-sigma ADCs (data rates 125 SPS – 8 kSPS)
- $8 \mu\text{V}$ pk-pk noise (PGA gain = 6, BW = 150 Hz)
- CMRR: -105 dB
- Integrated Right-Leg Drive amplifier
- Integrated Lead-off detection
- Integrated respiration impedance measurement (ADS1292R)
- Integrated test signals for verification
- Integrated low-drift ADC reference
- Integrated oscillator



Clinical ECG Portfolio



TIDA-010005/Software-Configurable Cardiac Pacemaker Detection Module Reference Design

Features

- Programmable Thresholds for detecting rise time, amplitude, duration and polarity of pace pulse
- Availability of measured parameters and conditioned pace signal on I²C Bus for further analysis and processing
- Interface with TI's ADS129X series of ECG front end devices
- Detection of negative pulse by handling polarity of pace pulse by software command

Applications


- Multi-parameter Patient Monitor
- Medical Sensor Patches
- Electrocardiogram (ECG)

Tools & Resources

- [TIDA-010005 and Tools Folder](#)
- [Design Guide](#)
- Design Files: [Schematics](#), [BOM](#), [Gerbers](#), [Software](#).
- Device Datasheets:
 - [INA317](#)
 - [TLV9062](#)
 - [TLV1702](#)
 - [DAC5578](#)
 - [ADS7142](#)

Benefits


- Detection of various types of pacemaker pulse provides the user one stop solution
- Modular approach compatible with other patient monitoring TI Designs such as TIDA-01614
- Compact form factor
- A flag signal and onboard LED are the indication of presence of 'valid' pace signal.





Design Challenges TIDA-010005 Solves

Design Challenge: Low noise, highly integrated signal path for ECG with lead-off and pace detection: to detect amplitude, rise time, width and polarity for diagnostics with different makes of Pacemaker

- Pace signals have specific characteristics:
 - Amplitude: 2 mV to 700 mV
 - Rise Time: 10 μ s to 200 μ s
 - Pulse Duration: 100 μ s to 2 ms
 - Polarity: +Ve or -Ve
- Major difficulty comes due to the small amplitude, narrow-width and varying slope of the pulse accompanied with noise consisting of: the ECG signal itself, noise from muscle movement due to beating of heart and breathing, 50/60 Hz pickup from the surroundings



Proposed Solution



In order to define a 'valid' pace pulse, the system discriminates various signals by comparing it to the user defined thresholds. The criteria for the same is:



- $\text{Amplitude} > V_{Amp_th}$
- $V_{RT_min} < \text{Rise Time Measure } (V_{RT}) < V_{RT_max}$
- $V_{PT_min} < \text{Pulse Duration Measure } (V_{PT}) < V_{PT_max}$

V_{Amp_th} , V_{RT_min} , V_{RT_max} , V_{PT_min} and V_{PT_max} comes from the DAC


TIDA-010005 Technical Details

Table 1. Key System Specifications

PARAMETER	SPECIFICATIONS	DETAILS
Input Voltage	3.3 V	Supply voltage to the system
PGAP	Analog non-inverting ECG input to system	Analog output of ADS129X PGA non-inverting pin which serves as input to the TIDA-010005
PGAN	Analog inverting ECG input to system	Analog output of ADS129X PGA inverting pin which serves as input to the TIDA-010005
Rise-time (T_R) measurement range	30–200 μ s	The measurement range of the rise time of the pace pulse
Pace pulse duration (T_D) measurement range	0.1–2 ms	The measurement range of the duration of the pace pulse
Input Pace signal amplitude range	8 mV – 700 mV	The range of amplitude of the pace signal tested with the TIDA-010005
Average steady-state current consumption	8 mA	Current consumption of the system
Clock A frequency	1.75 MHz	Onboard clock to the counter for T_R measurement
Clock B frequency	806 kHz	Onboard clock to the counter for T_D measurement
ADC resolution	12 bits (2 channel)	ADS7142 ADC to capture the rise time and duration
DAC resolution and channels	8 bits (8 channel)	DAC5578 sets the various threshold for pace pulse. I ² C enabled



Negative Pace Signal Detection - Negative Pace Pulse Undetected



Negative Pace Signal Detection Negative Pace Pulse Detected by setting the DACs Polarity Check (CH 6) Output High

Test Pacemaker Detection with TIDA-010005 & TIDA-01614



TIDA-01580 Wearable, Wireless, Multi-Parameter Patient Monitor Reference Design

Features


- Simple Wearable Multi-Parameter Patient Monitor for Photoplethysmography (PPG) and Electrocardiography (ECG)
- Provides Raw data to calculate heart-rate, Oxygen Concentration in Blood (SpO₂) and Pulse-transit Time (PTT)
- Uses Single-chip Bio-sensing Front-End AFE4900 for Synchronized ECG & PPG
 - PPG (Optical heart-rate monitoring and SpO₂) supports 4 LEDs and 3 PDs with Digital Ambient subtraction to improve the SNR
 - ECG (LEAD I) signals
- Integrated ARM Cortex-M3 + 2.4GHz RF Transceiver (CC2640R2F) supports wireless data transfer – BLE 4.2 and 5
- Operated from CR3032 (3V, 500mA Coin Cell Battery) with battery life of 30 days using highly efficient DC/DC converters
- Small form factor helps in easy adaptation to wearable applications

Target Applications

- Wireless Patient Monitor
- Wearable Fitness & Activity Monitor

Pulse Oximeter
ECG

Tools & Resources




Device Datasheets:

- [AFE4900 - CC2640R2F](#)
- [TPS61098 - TPS63036](#)
- [TPD1E10B06](#)

Benefits


- PPG supports 4 LEDs and 3 PDs with Digital Ambient subtraction to improve the SNR
- AC and DC lead off detection helps in correct measurement of vital signs
- Continuous Monitoring with lower operating power ensures battery life of 30 days
- Flexibility of ultra low power modes and integrated FIFO can keep MCU into sleep to increase the battery operation time




Copyright © 2018, Texas Instruments Incorporated

TIDA-01580 Medical Patches

This board is connected on bottom of the main board




This side is touching the wrist of one hand (ELECTRODE 1)



Other hand can be touched on PAD on top layer of the main board. (ELECTRODE 2)



Side View




- LAUNCHXL-CC2640R2F receives the signals remotely and displays on LabView GUI
- The design uses BLE 5.0 with an advertising time = 100ms

Design Challenges TIDA-01580 Solves

Design Challenge #1: Integration of multiple modalities at optimum SNR levels and small form factor



- Capturing synchronized ECG and PPG to enable PTT and BP calculations (non-invasive and without cuff)
- Pulse Transit Time (PTT): Time difference between the R-peak in the ECG waveform and the arrival of the blood pressure wave
- Simultaneous measurement of ECG and PPG together
- Along with other variables, such as the patient's size, weight, age, etc., algorithms show the correlation between PTT and systolic blood pressure.
- Challenging to synchronize both measurements – timing is the key!
(Powering up, clock timing, phase, drift with temperature)



Design Challenges TIDA-01580 Solves

Design Challenge # 2: BLE connectivity that does not interfere with measurement accuracy




- Signal amplitude: 0.2mV to 2mV (p-p)
- BW : as broad as 0.05 Hz to 300 Hz (Pace detection increases the bandwidth further)
- Reject environmental electrical signals, such as ac mains, security systems, and RFI to amplify and display the ECG signal
- Good CMRR of the signal chain and Right-leg drive (RLD) for CM rejection
- Differential- and common-mode filtering, environmental shielding, and algorithms



Design Challenges TIDA-01580 Solves


Design Challenge # 3: Extending battery life to multiple days to enable portability & wearability

- Selecting the right power architecture to enable extended battery life up to 24 hours (for rechargeable batteries) or 7 days (for primary cells)
- Powering with right buck, boost or buck-boost device instead of directly powering from battery (bypass modes in DC/DC converters)
- Sleep / shutdown / standby / deep sleep modes for radio devices like BLE, Wi-Fi etc.
- Selection of right battery charger (charging rate, termination current and quiescent currents are important!)



SHIP mode

- By definition: Ship mode electronically disconnects the battery from the rest of the system to minimize power drain while the product is idle. When the consumer turns the product on for the first time, the battery connects to the rest of the system and stays connected until the system decides to put itself back into ship mode.




SHIP mode



Out of SHIP mode

- Existing TI collaterals:
 - [Implementing Ship Mode Using the TPS22915B Load Switch](#)
 - [Don't let your battery drain on the shelf – use ship mode](#)
 - [Ship mode reference design \(TIDA-00556\)](#)



Low-voltage switchers

- For battery operated patches, a buck-boost and/or a boost is required. Few proposed devices are:

	TPS63802	TPS63030	TPS63036
Vin (Min) (V)	1.3	1.8	1.8
Vin (Max) (V)	5.5	5.5	5.5
Vout (Min) (V)	1.8	1.2	1.2
Vout (Max) (V)	5	5.5	5.5
Iout (Max) (A)	2	0.5	0.5
Rating	Catalog	Catalog	Catalog
Topology	Buck-Boost Buck/Boost	Buck-Boost	Buck-Boost
Type	Converter	Converter	Converter
Switching frequency (Min) (kHz)	2500	2200	2200
Switching frequency (Max) (kHz)	2500	2600	2600
Iq (Typ) (mA)	0.011	0.025	0.025
Features	Enable Light Load Efficiency Load Disconnect Power Good Pre-Bias Start-Up Synchronous Rectification UVLO Fixed	Enable Frequency Synchronization Light Load Efficiency Load Disconnect Synchronous Rectification UVLO Fixed	Enable Frequency Synchronization Light Load Efficiency Load Disconnect Synchronous Rectification UVLO Fixed
Switch current limit (Typ) (A)	4.5	1	1
Duty cycle (Max) (%)	100	100	100
Operating temperature range (C)	-40 to 125	-40 to 85	-40 to 85
Package Group	VSON-HR 10	SON 10	DSBGA 8

	TPS61022	TPS61099	TPS61291	LMR61428
Vin (Min) (V)	0.5	0.7	0.9	1.2
Vin (Max) (V)	5.5	5.5	5	14
Vout (Min) (V)	2.2	1.8	2.5	1.24
Vout (Max) (V)	5.5	5.5	3.3	14
Switch current limit (Typ) (A)	8	1	1	2.85
Regulated outputs (#)	1	1	1	1
Switching frequency (Min) (kHz)		1000	2300	300
Switching frequency (Max) (kHz)	1000	3000	2300	2000
Iq (Typ) (mA)	0.027	0.0008	0.005	0.1
Features	Enable Light Load Efficiency Load Disconnect Pre-Bias Start-Up Synchronous Rectification UVLO Fixed	Synchronous Rectification UVLO Fixed	Bypass Mode	Enable
Duty cycle (Max) (%)	97	90	90	70

Full system | Multi-parameter patient monitor + wireless sensors



Nurse's Station /
Doctor's Office

WiFi
CERTIFIED



Bluetooth®

Bluetooth®

Bluetooth®



Why TI SimpleLink™ for Multi-parameter patient monitor + sensor patch?

Low Power



- BLE SoC with integrated Ultra low Power Sensor Controller
- Wi-Fi low power IoT
- Best-in-class standby current

Ease of use




- CC3135/CC3235 Wi-Fi modules
- 5GHz Wi-Fi to reliably connect to hospital network
- BLE multi-role support
up to 32 simultaneous connections

Secure



- FIPS 140-Level 1 validation
- Offload CPU bandwidth – HW crypto accelerators
- Secure boot


Small size



- BAW: First crystal-less wireless BLE SoC
12% area savings in reference design
- Tiny BLE SoC: CC2640R2F
2.7mmx2.7mm DSBGA

Invest once, reuse effortlessly

100% code reuse



- Learn more about SimpleLink code portability
- SimpleLink Medical Resources

Bluetooth®

Full-featured Bluetooth 5 solutions

Features all mode Bluetooth 5 certified support, automotive-qualified wireless MCUs, and has the industry's smallest full-featured Bluetooth 5 solution.

[BLE overview](#)

[View our BLE products](#)

Sub-1GHz

Ultra-low power long-range star network

Out-of-box star network solution, multi-year operation on coin cell battery. Programmable ultra-low power sensor controller interface.

[Sub-1GHz overview](#)

[View our Sub-1GHz product](#)

WiFi CERTIFIED

Dual-core and low power SoCs

FIPS-verified ICs optimized for low power. Enhanced application and network security capabilities.

[Wi-Fi overview](#)

[View our Wi-Fi products](#)

zigbee

Certified mesh network

Zigbee 3.0-certified with lowest power integrated +20dBm PA and green power support.

[Zigbee overview](#)

[View our Zigbee products](#)

THREAD

Self-healing low-power mesh network

Lowest power thread platform. OpenThread stack support. Optimized router examples available in SimpleLink academy.

[Thread overview](#)

[View our Thread products](#)

((o)) Multi-standard

Concurrent wireless protocol operation

Concurrent multi-protocol & multi-band. BLE + Zigbee or BLE + Sub-1 GHz. Pre-built multi-protocol manager with flexible priority scheduler.

[Multi-standard overview](#)

[Multi-standard products](#)

Common software

SDK

TIDA-010029

Wearable, 16-phase multi-sensor SpO₂ and heart rate monitor (HRM) reference design with Bluetooth® 5

Features

- Provides raw data to calculate heart rate, SpO₂, and other related parameters
- Uses single-chip, bio-sensing, front-end AFE4420 device for PPG measurement
 - Enables signal acquisition of up to 16 phases and multi-wavelength measurements with the flexible allocation of LEDs and photodiodes in each phase
- Integrated Arm® Cortex®-M3 and 2.4-GHz RF transceiver (CC2640R2F) supports wireless data transfer through Bluetooth® low energy 4.2 and 5.0

Applications

- Pulse Oximeter
- Wearable Fitness and Activity Monitor
- Multiparameter Patient Monitor
- Medical Sensor Patches


Tools & Resources

- TIDA-010029 and/or Tools Folder**
- Design Guide**
- Design Files:** Schematics, BOM, Gerbers, Software, etc.

- Device Datasheets:**
 - AFE4420
 - CC2640R2F
 - TPS63036
 - TPS61099
 - TPD1E10B06

Benefits

- PPG supports 4 LEDs and 4 PDs with Automatic Ambient subtraction to improve the SNR
- Continuous Monitoring with lower operating power ensures battery life of 100 hours for continuous operation
- Flexibility of ultra low power modes and integrated FIFO can keep MCU into sleep to increase the battery operation time
- Small form factor helps in easy adaptation to wearable applications




Copyright © 2018, Texas Instruments Incorporated

Design Challenges TIDA-010029 Solves

Design Challenge # 1: Integration of multiple modalities at optimum SNR levels and small form factor

- Provides raw data to calculate heart rate, SpO₂, and other related parameters
- Supports 4 LEDs and 4 PDs with improved the SNR
- External Light sources can disturb the receive signal (Smart automatic ambient light suppressing technologies – integrated into the device)
- Standard pulse oximetry uses Red and IR lights. The additional green and orange wavelengths can distinguish to the absorption of deoxyhemoglobin (RHb) and oxyhemoglobin (HbO₂).
- Physiology of the Human body (Skin tone and body Location determines use case, Transmissive vs reflective oximetry)






Source: An applicable approach for extracting human heart rate and oxygen saturation during physical movements using a multi-wavelength illumination optoelectronic sensor system ([Link](#))




Design Challenges TIDA-010029 Solves

Design Challenge # 2: BLE connectivity that does not interfere with measurement accuracy

- Supports wireless data transfer through Bluetooth® low energy 4.2 and 5.0 with CC2642R
- CM chokes placed at the PD inputs to reduce noise effects on very sensitive signals coming from photo diodes
- Layout Optimization for AFE and BLE to minimize interface while maintain the measurement accuracy
- AFE4420 and CC2640R2F placed on different layers with GND layer in between for shielding
- Split ground for RF and AFE sections
- STAR grounding at the negative terminal of the battery



TIDA-010029 board details



TIDA-010043: Efficient, high-current, linear LED driver reference design for SpO₂ and other medical applications

Features

- Higher efficiency by providing headroom control to optimally drive LED
- Diagnostics (LED open, LED short, LED disconnect) and protection mechanism for fault by disconnecting the supply using e-fuse
- Wide operating input range from 1.8 V to 5.5 V to support all battery type inputs
- Programmable range selection for LED current (from 100 mA to 1.5 A)
- Compatible with AFE44xx series EVMs (tested with AFE4403 EVM) for driving LEDs
- Tested with MSP430G2553 LAUNCHPAD™ and two on-board LEDs (green and red)

Target Applications

Pulse oximeter
monitor

Multiparameter patient

Tools & Resources




Device datasheets:

- [TPS63802](#)
- [CSD13306W](#)
- [TMUX1119](#)
- [OPA333](#)
- [OPA365](#)
- [TPS2595](#)
- [CSD1320Q2](#)
- [TLV7021](#)

Benefits





- Built-in self-calibration and look-up table for headroom control
- Dynamic adjustment of the LED supply voltage, to achieve optimal efficiency at the programmed current levels




Design Challenges TIDA-010043 solves

Design Challenge: Accurate SpO₂ measurement across varying photodiode placements, wavelengths and body physiologies through optimizing LED driving circuit and photo-diode signal chain at low power

- External Light sources are disturbing the receive signal (Smart ambient light suppressing technologies – integrated into the device)
- Physiology of the Human body (Skin tone and body Location determines use case, Transmissive vs reflective oximetry)
- Optical Signal Path (Define optimal Photo Diode size & LED wavelengths, use multiple sensors and PDs, optimum in PD / LED spacing)




TIDA-010043 Technical Details – flowchart




TIDA-010043 Technical Details – Calibration

Software steps for calibration (the current through LED will increase by 100mA for each step)



TIDA-010043 Technical Details – block diagram



Calibration loop

VCC = 5 V (from LaunchPAD™ or AFE EVM)

TIDA-01624 Bluetooth-Enabled High Accuracy Skin Temperature Measurement Flex PCB Patch

Features

- High Accuracy, Low Power Temperature Sensor
- BLE 4.2 and 5 enabled microcontroller
- Thin-Film Flexible Battery Power, enabling entirely flexible design
- Integrated PCB antenna
- Temperature updates every second

Benefits


- Low power consumption and long battery life
- Extremely long shelf life (3+ Years)
- Small, Flexible Form Factor
- Connects to Smart Device
- Zero-Calibration to $\pm 0.1^{\circ}\text{C}$ Accuracy

Applications

- Medical Sensor Patches
- Multiparameter Patient Monitors
- Smart Patches

Tools & Resources

- TIDA-01624 and/or Tools Folder
- Design Guide
- Design Files: Schematics, BOM, Gerbers, Software, etc.
- Device Datasheets:
 - TMP117
 - CC2640R2F



TMP117x

Ultra-High Accuracy Digital Temp Sensor with integrated non-volatile memory

Features

Accuracy


TI Part	Accuracy (°C)	Accuracy Full Range
TMP117M	±0.1°C @ (30°C to 45°C)	±0.2°C @ (0°C to 85°C)
TMP117	±0.1°C @ (-20°C to 50°C)	±0.3°C @ (-55°C to 150°C)
TMP117N	±0.2°C @ (-40°C to 100°C)	±0.3°C @ (-55°C to 150°C)

- 16-bit Resolution (0.0078°C)
- Minimum PSRR: 1LSB = 7.8 m°C/V

Integrated EEPROM

Low power consumption

- 140µA Iq during conversion
- 3.5µA Average Iq @ 1Hz
- 150nA Shutdown Iq
- 1.8V – 5.5V



Digital feature: Automatic offset NVM/ Soft Reset

Interface: Single wire

Packaging

- 6pin WSON (2 x 2) mm
- 6pin WCSP (1.6 x 1) mm

Applications

- Gas Meter
- Medical
- Cold Chain
- Wearables
- Instrumentation & Test
- Thermocouple – Reference

Benefits

Ultra-high Accuracy

- Meets ASTM E1112 & ISO medical standards:
 - 0.1°C acc. range 35.8°C to 42°C
 - No calibration needed; NIST Traceable



Integrated Non-volatile memory

- Store configuration even after losing power
- 64 Bits of general purpose scratch pad memory



Low Power Consumption

- 3.5µA Average Iq @ 1Hz; serial bus inactive
- 150nA Shutdown Iq; serial bus inactive

Digital feature & I2C Interface:

- Programmable Temperature Alert & Offset value
- Soft Device Rest

Smallest Package: 6 PIN, QFN & CSP



TIDA-010040 / Alarm tone generator reference design

Features

- Provides a solution that can be used for producing auditory alarms described in the IEC60601-1-8 medical specification.
- Capable of outputting low, medium, and high priority alarms with software-adjustable rise/fall time, pulse duration, pulse spacing, and burst spacing
- Can output 8 different alarm melodies as described in the IEC60601-1-8 spec.
- Coincidence Detector confirms if speaker is actually making correct sound


Benefits

- Timing parameters are all adjustable via firmware
- Minimal code-space used allows for low-cost MSP430
- Booster pack design makes for easy setup

Target Applications

- [Multiparameter Patient Monitor](#)
- [Dialysis Machine](#)

[Infusion Pump](#)
[Surgery Equipment](#)







Key Design Challenge: Technical details

IEC60601-1-8 based medical alarm tone generation & coincidence detection



- Provides a solution that can be used for producing auditory alarms described in the IEC60601-1-8 medical specification.
- With a Low cost microcontroller and some external hardware, Analog Burst patterns are created. The Rise and Fall times, the Amplitude ,width and Frequency of these Bursts are all programmable.
- The Alarm made by the circuit when attached to a Patient monitor might be indicative of a Medical condition such as a Heart attack. Therefore it is very important that the Alarm itself has some kind of monitoring of its operation.

A microphone is Placed close to the speaker and it monitors the sounds made by the Alarm. If the Actual sound received by the microphone is not the same as the sound intended by the circuit an error flag is raised. This is called Coincidence detection



Coincidence
Detection

TIDA-010040 – Coincidence Detection



Welcome to the TI E2E™ Community

Find answers to your technical questions

 Search E2E by part number and/or keyword. (e.g. OPA333 output peaking)

2,785 Contributing TI employees 295,218 Issues resolved

Recent activity

Forums	Blogs
--------	-------

TLV320AIC34EVM-K: Audio Serial Data Bus

Part Number: TLV320AIC34EVM-K Hi all, I have to read two stereo channels, so two I2S outputs A and B of the Codec are used. Both I2S_A and B outputs are connected to an FPGA on two IP I2S. Is it possible to read I2S_A (Audio Serial Data Bus A) and I2S_B (Audio Serial Data Bus B) outputs simultaneously ? That means both I2S outputs are independent, i.e no data multiplexed ? Best regards, Pa...



Audio forum

Patopat23

10/22/2018 6:37:32 AM

0 Views | 0 replies

Top contributors



YiKai Chen
574855 points



Jens-Michael Gross
227195 points



IMPORTANT NOTICE AND DISCLAIMER

TI PROVIDES TECHNICAL AND RELIABILITY DATA (INCLUDING DATASHEETS), DESIGN RESOURCES (INCLUDING REFERENCE DESIGNS), APPLICATION OR OTHER DESIGN ADVICE, WEB TOOLS, SAFETY INFORMATION, AND OTHER RESOURCES "AS IS" AND WITH ALL FAULTS, AND DISCLAIMS ALL WARRANTIES, EXPRESS AND IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT OF THIRD PARTY INTELLECTUAL PROPERTY RIGHTS.

These resources are intended for skilled developers designing with TI products. You are solely responsible for (1) selecting the appropriate TI products for your application, (2) designing, validating and testing your application, and (3) ensuring your application meets applicable standards, and any other safety, security, or other requirements. These resources are subject to change without notice. TI grants you permission to use these resources only for development of an application that uses the TI products described in the resource. Other reproduction and display of these resources is prohibited. No license is granted to any other TI intellectual property right or to any third party intellectual property right. TI disclaims responsibility for, and you will fully indemnify TI and its representatives against, any claims, damages, costs, losses, and liabilities arising out of your use of these resources.

TI's products are provided subject to TI's Terms of Sale (www.ti.com/legal/termsofsale.html) or other applicable terms available either on ti.com or provided in conjunction with such TI products. TI's provision of these resources does not expand or otherwise alter TI's applicable warranties or warranty disclaimers for TI products.

Mailing Address: Texas Instruments, Post Office Box 655303, Dallas, Texas 75265
Copyright © 2019, Texas Instruments Incorporated